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ANSI ASHRAE

## SUMMARY

<b>1</b>	<b>GENERAL INFORMATION .....</b>	<b>4</b>
1.1	TECHNICAL SUPPORT.....	4
1.2	USER TRAINING .....	4
1.3	GRAPHIC SYMBOLS.....	5
1.4	WARNING .....	6
1.5	RISK ASSESSMENT .....	7
1.6	STANDARD AND CERTIFICATION .....	8
<b>2</b>	<b>TRANSPORT AND STORAGE .....</b>	<b>9</b>
2.1	TRANSPORT CAUTION .....	9
2.2	STORAGE CAUTIONS .....	9
<b>3</b>	<b>PRODUCT DESCRIPTION .....</b>	<b>10</b>
3.1	PURPOSE.....	10
3.2	TYPICAL WORKING.....	10
3.3	GENERAL VIEW .....	11
<b>4</b>	<b>TECHNICAL DATA.....</b>	<b>12</b>
4.1	DIMENSIONS AND WEIGHTS .....	12
4.2	ELECTRICAL DATA.....	13
4.3	WORKING CONDITIONS .....	13
<b>5</b>	<b>INSTALLATION .....</b>	<b>14</b>
5.1	USEFUL SPACE AND CROSSING .....	14
5.2	ROOM AIR MAKE-UP .....	15
5.3	INSTALLATION.....	16
5.4	HYDRAULIC FITTING / GASES CONNECTION.....	17
5.5	ELECTRIC CONNECTIONS.....	18
5.5.1	ELECTRIC PANELS ON FUME HOOD .....	19
5.5.2	MOTOR FAN .....	20
5.6	EXHAUST FUMES DUCTING .....	21
5.7	IN-LINE FILTERING.....	21
5.8	AUTOMATIC COMPENSATING DEVICES / MULTIPLE FUME HOODS.....	21
5.9	FALSE CEILING AND DIFFUSERS .....	22
5.10	WARNING AND PERIODIC CONTROL .....	22
<b>6</b>	<b>OPERATING .....</b>	<b>23</b>
6.1	START-UP .....	23
6.2	VARIABLE AIR VOLUME SUCTION CONTROL .....	23
6.2.1	V.A.V. SUPPLY AND INSTRUMENT LINKS.....	24
6.2.2	V.A.V. APPARATUS.....	25
6.2.3	OPERATING.....	25
6.2.4	NORMAL MODE.....	26
6.2.5	BASE PARAMETER.....	26
6.2.6	WORKING PARAMETERS .....	27
6.2.7	WORKING PARAMETERS SETTING.....	28
6.2.8	ALARM MANAGEMENT.....	29
6.2.9	OPTIONAL PERFORMANCES .....	30
6.3	INSTRUMENT CALIBRATION.....	31
6.3.1	CALIBRATION MENU .....	31
6.3.2	“CAL” - FRONTAL AIR VELOCITY MEASUREMENT CALIBRATION .....	32

6.3.1	“OFF” – AIR VELOCITY MEASUREMENT COMPENSATION .....	33
6.3.2	“UEL” – INVERTER MINIMUM FREQUENCY RATE.....	33
6.4	AIR VELOCITY MEASURING PROBE .....	34
6.5	INVERTER SSD DRIVES - 650 SERIES.....	34
6.6	CIRCULAR SHUT-OFF DAMPERS – ELECTRONIC ACTUATOR .....	35
6.7	FUMES SUCTION AND EXHAUST .....	35
6.8	SASH.....	36
<b>7</b>	<b>CARE AND MAINTENANCE .....</b>	<b>37</b>
7.1	GENERIC ADVICES .....	37
7.2	WORKING AREA ACCESS .....	37
7.3	SASH AND COUNTERWEIGHT MAINTENANCE .....	37
7.4	TO AVOID .....	38
<b>8</b>	<b>SPARE PARTS .....</b>	<b>39</b>
<b>9</b>	<b>ELECTRIC diagrams.....</b>	<b>40</b>
9.1	AIR SUCTION CONTROL DIAGRAM – POWER UP TO 0,75 kW.....	40
9.2	AIR SUCTION CONTROL DIAGRAM – POWER UP TO 1,5 kW.....	41
9.3	GENERAL BLOCK DIAGRAM .....	42
<b>10</b>	<b>SUPPLYING FIXTURE .....</b>	<b>43</b>
<b>11</b>	<b>CERTIFICATIONS.....</b>	<b>45</b>
11.1	“CE MARK”.....	45
11.2	EN 14175 - PART 1,2,3 CERTIFICATION.....	46
11.3	EN 14175 - PARTE 6 (V.A.V.) CERTIFICATION.....	47
11.4	ANSI/ASHRAE 110-1995 CERTIFICATIONS.....	48

## 1 GENERAL INFORMATION

### 1.1 TECHNICAL SUPPORT

---

**MANUFACTURER :****BICASA Srl**

Viale delle Industrie, 33  
ITALY – 20044 Bernareggio (MB)  
Ph. +39 039 60291  
Fax +39 039 6093153  
[www.bicasa.it](http://www.bicasa.it) - [info@bicasa.it](mailto:info@bicasa.it)

**MAIN OFFICE & SERVICE :****BICASA Srl**

Viale delle Industrie, 33  
ITALY – 20044 Bernareggio (MB)  
Ph. +39 039 60291  
Fax +39 039 6093153  
[www.bicasa.it](http://www.bicasa.it) - [info@bicasa.it](mailto:info@bicasa.it)

**FRENCH OFFICE :****BICASA FRANCE Sarl**

3, rue des Artisans  
France - 68600 Algsheim  
Ph. +33 3 89727063  
Fax +33 3 89221756  
[bicasafrance@bicasa.eu](mailto:bicasafrance@bicasa.eu)

### 1.2 USER TRAINING

---

The information and data content in the present manual are to be generically considered and without any scientific effect. Safety laws are particularly rigorous and restrictive; we suggest therefore a careful preliminary estimation.

This apparatus has been designed to CHEMICAL multipurpose applies: performances are unidentified. This fume hood is not able to BIOLOGY an BIO-HAZARD applications.

Manufacturer's responsibility is limited to the original installed components, since it is not possible to check the effective use and the possible modification which would be carried out after delivery. Every change (even not structural or operative) release the manufacturer from the responsibility of eventual failures, damages, malfunctions etc., further to generate the immediate decadence of warranty, maintenance and certification.

This manual must be kept available for operator's consultation and for the staff entrusted for fume hood setting up. Installation and ordinary maintenance are to be performed by manufacturer's personnel or authorized technicians. During installation, hydraulic and electric plants are to be set up by skilled and qualified technicians, in compliance with the regulation set forth. The technician would release a Conformity Certificate of the installed hydraulic/electric plant. The use of the fume hood is restricted to trained personnel who should however carefully read the present manual. The fume hood is solely destined to protect the operator from fumes or chemical vapours; it is not therefore suitable for protecting against pathogen or radioactive compounds. Before starting any work, the operator should estimate if the equipment complies with his specific requirement. The fume hood is not included in the medical devices category and, because its use is restricted to the laboratory and because there are not moving parts operating by electrical or mechanical device, it is not subjected to CE machines mark.

Should any doubts or misunderstandings regarding the correct use of the unit be arise, contact the manufacturer or the service centre nearest to you.

Upon request, the “AFTER INSTALLATION TEST” service may be proposed, according to EN 14175 – part 4 standard.

In order to get any other necessary information related to this equipment, please contact your supplier or the manufacturer directly to <http://www.bicasa.it>:

### 1.3 GRAPHIC SYMBOLS

---

**GENERIC DANGER** (look at manual instruction)

---

**ELECTRIC DANGER** (look at manual instruction)

---

**FIRE DANGER** (look at manual instruction)

---


**CORROSIVE SUBSTANCES** (look at manual instruction)

---

**BE CAREFUL** (look at properly working mode)

---

**NOTE** (helpful information and recommending)

---

**PROTECTION HEART CONNECTION**

---

**EQUIPOTENTIAL CLAMP**

---

**SPECIAL SETTING SCRAP**

---

## 1.4 WARNING

---

Alteration, manipulation or modification of single component must **ALWAYS** be **AVOIDED**, even if considered of limited importance. The ordinary maintenance and part replacement are to be carried out by manufacturer's personnel or skilled technicians authorized by the manufacturer.

The manufacturer declines any responsibilities for damages occurred to people or things due to equipment misuse.

Access to tools and electrical components (even if low-voltage operating) is prevented by suitable mechanical safety devices (closing device which require specific mechanical tool). Any operation practiced on the electrical plant, must be carried out by manufacturer's personnel or skilled technicians authorized by the manufacturer

Before approaching any kind of electrical equipment, make sure that power supply shut off. Before approaching any kind of fittings, particularly combustible gas lines, make sure that main feed valves was **CLOSED**.

A careless use, the access to the electrical components, the manipulation of gas and fluid feeding lines, the removal of external coverage may cause **HIGH DANGER RISK**.

This kind of fume hood is suitable for **CHEMICAL** use: it cannot therefore be used as protecting instrument against biological hazards.

Should any doubt or misunderstanding regarding the correct use of the unit be arise, contact the manufacturer.

It is very **DANGEROUS** the unwary or imprudent use of the hood; the covering panel removal and the access to inside components; the alteration of fixtures, fittings and piping.

## 1.5 RISK ASSESSMENT

---

Fume cupboard is designed for chemical hazardous application. The multipurpose apply does not allow to manufacturer to estimate final process risk, even when he knows generically the final destination.

**Operator healthcare and safety are always owner responsibility.**

Estimating the process and the apply risks, the operator has to consider:

- Kind of activity;
- The material under analysis or manipulation;
- The potential substance alteration and reaction;
- Processing-time and hazard exposure;
- Residue remained inside the working chamber;
- Working area environmental factor influence.

Before starting the fume hood, must be taken care to the latent hazard inference, properly upgrading procedure and time test in compliance with the effective equipment protection capacity. Try to reduce as much as possible the risk by using simple shrewdness like:

- keep sash always as closed as possible (when the fume hood is not equipped with automatic sash closing device);
- reduce essentially reagent and sample quantities;
- check carefully process temperatures, limiting them as much as possible;
- reduce fumes and gas development by using cold-condensers when possible;
- carefully estimate potential fire or explosion risk, handling flammable or explosive compounds and reagents.

**IT IS RIGOROUSLY FORBIDDEN TO USE THE HOOD WITHOUT VENTILATION,  
THUS WITHOUT MECHANICAL CONVECTION.**

**The robustness performance of the hood will stop without forced ventilation.**

It is essential to verify the correct air flow efficiency on working volume referred to performance carried out:

- verify incoming air flow rate, measuring the frontal air velocity by an anemometer gauge (opening sash at cm 45 of height);
- calculate the working volume number of change.


## 1.6 STANDARD AND CERTIFICATION

LITE fume cupboards are in compliance with the following standards:

- **EN 14175 : 2004 Part 1, 2 and 3**, European standard – Fume cupboards.  
“Vocabulary; Safety and performance requirements; Type test methods”
- **EN-14175 : 2006 Part 6**, European standard – Fume cupboards.  
“Variable air volume fume cupboards”.

- **ANSI-ASHRAE 110-1995**, American standard  
“Guideline, Method of Testing Performance of Laboratory Fume Hood”.


LITE BT 1.2, LITE BT 1.5 and LITE BT 1.8 had been tested by independent and qualified bureau of certification. Certification sheet copy attached to this manual.

LITE fume cupboard equipment and apparatus according to:

<b>EN 61010-1 : 2001 (2<sup>nd</sup> Edition)</b>	Safety requirements for electrical equipment for measurement, control and laboratory use - Part 1: General requirements
<b>EN 61326-1/EC</b>	Electrical equipment for measurement, control and laboratory use – EMC requirements – Part 1: General requirements
<b>2004/108/CEE</b>	Electromagnetic compatibility
<b>2006/95/CE</b>	Low voltage
<b>EN 60529</b>	Degrees of protection provided by enclosures (IP code)
<b>CEI 23-50</b>	Plugs and sockets for domestic use
<b>EN 60309-1</b>	Plugs and sockets for industrial use - Part 1
<b>EN 60309-2</b>	Plugs and sockets for industrial use - Part 2
<b>EN 60947-2</b>	Low voltage switch gear and control gear. Circuit-breakers
<b>EN 60898</b>	Circuit-breakers for over current protection for household and similar installations
<b>FG7OM1 LS0H</b>	Low Smoke Zero Halogen labeled electric cables
<b>EN 12150-1</b>	Safety glasses according to E3 European standard
<b>EN 13792</b>	Color coding of taps and valves for use in laboratories
<b>DIN 12918</b>	Taps and valves for use in laboratories

Electrical fans, if fitted with fume cupboard, are in compliance with 98/37/CEE:1998 and updating.

See also the attached CE declaration of conformity.

## 2 TRANSPORT AND STORAGE

### 2.1 TRANSPORT CAUTION

The LITE chemical fume hood is dismantled normally shipped, in order to reduce delivery volume and weight. Before being shipped, single components are adequately packed and protected for road, rail or boat container transport. Specific protection packages are provided for special or sea freight shipment.

In case of customer transport carry out, please get particularly care to glass parts, metallic polished surfaces, counter top and sinks, holding-out edge. It is appropriate to apply to specialized packing company.

Previous to remove a built-in hood, disconnect all supply net fittings (electrical, gases, water), disassemble the main frame components, protect them with adequate material, specially the most fragile parts as:

<b>GLASSES</b>	Strong package (wooden crate, rigid cardboard). Special care must be taken to the edges
<b>PAINTED ELEMENTS</b>	Wrap with polythene film to preserve the parts from the moisture. External shock-proof packing with light cardboard.
<b>PLASTIC LAMINATES</b>	Light cardboard package, with special care for the edges.
<b>FITTINGS, ELECTRIC PANELS, ACCESSORIES</b>	Wrap with polythene film to preserve the parts from the moisture. External shock-proof packing with light cardboard.

### 2.2 STORAGE CAUTIONS

Care must be taken to store the LITE fume hood or its components in covered and protected warehouses, with 75% RH max and temperatures between + 10°C to + 35°C.

When the fume cupboard taking out from store and before getting new installation, please inspect carefully the equipment to ensure that are undamaged, particularly electric and gas tubing conditions, sash and counterweight ropes integrity.

**DO NOT LEAVE the fume hood exposed to weather conditions**

Should any doubt or misunderstanding concerning transport be arise, contact the manufacturer or the service centre nearest to you.

## 3 PRODUCT DESCRIPTION

### 3.1 PURPOSE

The LITE fume hood is equipment tailored for fumes, toxic gases and hazardous vapors containment, further to protecting the operator and working room environment from chemical hazard.

The apparatus includes an aspiration working chamber with possibility to make any kind of chemical hazard process in presence of air forced suction.

Air flow is usually obtained by independent and remote suction fan (not on-board) separately purchasable in order to lay-out or by centralized air suction plant.

**This type of fume hood is EXCLUSIVELY destined to protect the operator from chemicals; it is not therefore suitable for protecting against biological or nuclear lab hazard.**

**LITE fume cupboard IS NOT AVAILABLE for NUCLEAR or RADIOACTIVE working, nor for BIOLOGICAL or BIOHAZARD (Laminar Flow cabinet) manipulation.**

### 3.2 TYPICAL WORKING

It is essential to carry out a preliminary project valuation about air flow system and plant process parameter. Without an appropriate performance of suction fan and/or air connecting pipes, the fume hood correct operating will be defaulting. A wrong air extraction process and connection will generate bad performance and high risk for operator safety.

The forced suction air flow performed by the remote fan or by central air suction system comes through the frontal opening, between sash-handle and countertop air-foil.

The air flow moves across inside chamber volume, channeled to special air collecting device up to exhaust duct on the hood ceiling. Afterward, vapors and fumes are ejected outside, using as well as possible molecular or carbon filtering case, separately purchasable in accordance with plant lay-out.

**Hood must be ALWAYS provided with an air exhaust system toward the outside. Although the suction piping equipped with in-line filtering apparatus, DO NOT RECYCLE air into the room.**



The aerodynamics shape of the sash handle and the air-foil presence on the countertop optimize the air flow cross on the frontal opening, ensured operating chamber volume exchange robustness.

**To assure the maximum level of safety for the period of process, the sash MUST TAKE-DOWN as soon as possible.**

### 3.3 GENERAL VIEW

## 4 TECHNICAL DATA

### 4.1 DIMENSIONS AND WEIGHTS

LITE fume cupboard made for each product version according to EN 141750 European Standard:

- **BENCH TOP :** with countertop at cm 90 of height and SINGLE sash, with sliding glasses
- **LOW LEVEL :** with countertop at cm 50 of height and DOUBLE sash with sliding glasses
- **WALK-IN :** without countertop (floor standing) and TRIPLE sash with sliding glasses

LITE - BENCH TOP			LITE BT 1.2	LITE BT 1.5	LITE BT 1.8	LITE BT 2.4
Outside Dimensions	(L/P/H)	mm	1200/900/2480	1500/900/2480	1800/900/2480	2400/900/2480
Inside Dimensions	(L/P/H)	mm	1140/650/1380	1440/650/1380	1740/650/1380	2340/650/1380
Height (max)	(H)	mm	2800	2800	2800	2800
Countertop (useful)	(L/P)	mm	1140/720	1440/720	1740/720	2340/720
Sash Opening (useful)	(H)	mm	900	900	900	900
Exhaust Duct	(Ø)	mm	200	250	250	315
Weight (*)		kg	260	300	360	460

(\*) Including gres countertop.

LITE - LOW LEVEL			LITE LL 1.2	LITE LL 1.5	LITE LL 1.8	LITE LL 2.4
Outside Dimensions	(L/P/H)	mm	1200/900/2480	1500/900/2480	1800/900/2480	2400/900/2480
Inside Dimensions	(L/P/H)	mm	1140/650/1780	1440/650/1780	1740/650/1780	2340/650/1780
Height (max)	(H)	mm	2800	2800	2800	2800
Countertop (useful)	(L/P)	mm	1140/720	1440/720	1740/720	2340/720
Sash Opening (useful)	(H)	mm	1400	1400	1400	1400
Exhaust Duct	(Ø)	mm	200	250	250	315
Weight (*)		kg	270	310	375	490

(\*) Including gres countertop.

LITE - WALK-IN			LITE WI 1.2	LITE WI 1.5	LITE WI 1.8	LITE WI 2.4
Outside Dimensions	(L/P/H)	mm	1200/900/2480	1500/900/2480	1800/900/2480	2400/900/2480
Inside Dimensions	(L/P/H)	mm	1140/650/2280	1440/650/2280	1740/650/2280	2340/650/2280
Height (max)	(H)	mm	2800	2800	2800	2800
Sash Opening (useful)	(H)	mm	1800	1800	1800	1800
Exhaust Duct	(Ø)	mm	250	250	315	350
Weight		kg	270	310	375	490

## 4.2 ELECTRICAL DATA

SUPPLY :	3/N/T 400V AC - 1/N/T 230V AC
FREQUENCY :	50/60 Hz
MAX POWER CONSUMPTION :	6,0 kW (3/N/T 400V AC) 3,5 kW (1/N/T 230V AC)
DEGREE OF PROTECTION :	IP55
GENERAL APPARATUS BREAKER : (frontal position – breaking capacity 6KA, type C)	Automatic switch 16A
SINGLE SOCKET MAX POWER CONSUMPTION : (interlocked socket – on request – with different input power marked on)	1000 W
LIGHTING :	≥ 1000 Lux

### EQUIPOTENTIAL EARTH CONNECTION

The hood frame must be **MANDATORY connected** to the equipotential earth net, using a adequate cable (4 mm<sup>2</sup> minimum) to connect to terminal junction box placed under the countertop (see 5.5 subsection “5.5 - ELECTRIC CONNECTIONS”).

### EXHAUST FAN

For correct use, the hood frame must be **MANDATORY channeled** to a remote fan or to a centralized air suction plant. Processing features would change depending on suction plant layout configuration and operating flow rate request on frontal hood opening.

See and ask for motor fan and ventilator documentation.







## 4.3 WORKING CONDITIONS

OPERATING TEMPERATURE	from +5°C to +35°C
OPERATING HUMIDITY (MAX)	75%

## 5 INSTALLATION

### 5.1 USEFUL SPACE AND CROSSING

Fume cupboard safety and correct operating needs an appropriate site allocation. According to CEN/TS 14175 standards Part 5 sub-chapter 4.2 with the exception of different National regulations indicate otherwise, the following notes recommend spaces and crossing ways between fume hood and laboratory elements or equipment. All dimensions given are for guidance to operators safety and protection as well as working efficiency

<p><b>MINIMUM DISTANCE HOOD – OPERATING ZONE:</b></p> <p>Minimum distance from hood sash to any part of laboratory frequently used by other personnel (excepted for wall)</p> 	<p><b>MINIMUM DISTANCE HOOD – BENCH :</b></p> <p>Minimum distance from hood sash to bench top opposite to it</p> 	<p><b>MINIMUM DISTANCE HOOD – WALL :</b></p> <p>Minimum distance from hood sash to a structured wall</p> 
<p><b>MINIMUM DISTANCE HOOD – DOORWAY :</b></p> <p>Minimum distance from hood sash to a doorway (excepted fire emergency exit)</p> 	<p><b>MINIMUM DISTANCE HOOD – EMERGENCY DOOR:</b></p> <p>Minimum distance from hood sash to a doorway (excepted fire emergency exit)</p> 	<p><b>AIR FLOW INTERFERENCE</b></p> <p>Carefully consider the air flow interaction between face to face sash hood</p> 

A ceiling height of 3 meters is recommended and height of 2.7 meters should be the minimum feature.

With lowered ceiling, particular hood models are available.

The possibility of fire or explosion in a fume cupboard should always be considered when positioning it into laboratory. Special attention should be given to the location of equipment and to national regulation about MANDATORY escape route from laboratory, doorway and emergency door..

## 5.2 ROOM AIR MAKE-UP

---

The fume hood operating causes an air consumption. The ambient where the fume hood is installed should therefore supply an air volume equivalent to the exhausted one (air make-up).

It is therefore essential assure an adequate air-inlet in the laboratory room, in order to respect both process and exhaust parameter. An air make-up (inside large room such as open space) may be happened naturally or by means of dedicated compensation grids (pass-through wall or pass-through door) or by forced air inlet throughout an oversized laboratory conditioning plant. The provisional fume hood operation must be however considered.

**Air shortage within the room means a decrease of fume hood performances and less safety for the operator.**

When laboratory activity requires and building condition allows the best way is to restore or to increase exhausted is by the HVAC system modification.

In large room the offset between exhaust air and make-up air is quit not necessary (e.g. open space lab).

In small room it could be restored by gratings installed on doors or on wall, linking other rooms or hallway.

Daily and weekly time-frame use of fume cupboard will be evaluated.

**Fume cupboard electronic master board equipped with free contact to interface the hood control system to HVAC.**

Large isolated obstructions at the side of the fume cupboard and projecting beyond the plane of the sash, e.g. an architectural column, can influence the performance of the equipment, particularly with V.A.V. system hood supplied.

The impact of room air inlet (location and dimension of) on the performance of fume cupboards should carefully considered: it could impair the regular working.

**According to EN-14175-4 - 5.8.3 standard, normally room air flow should not exceed 0.2 m/s at a zone 400 mm from the sash position.**

### 5.3 INSTALLATION

---

Fume cupboards are disassembled normally shipped. Installation process will start on site with fittings and fixtures designed.

When fume hood or components would be accumulated in temporary deposit prior to installation, check storage space and environment condition in order to minimize the risk of damage or deterioration. Care should be taken in protecting material stored before and during decoration and fitting out of the laboratory.

The adequacy of access to the laboratory should be inspected.

In difficult case, provisional elevator or scaffold access through windows may be necessary. If it would be, please consult specialized moving company supplied with adequate equipment and qualified staff.

Particular attention should be paid to:

- Slip road typology and dimension;
- Stairway and landing;
- Lift and freight elevator load capacity and dimension;
- Door dimensions;
- Corridor dimensions, change of corridor direction and potential obstruction.

**LITE fume hood should not be installed in laboratory where the environment condition would be as different as max RH of 75% and temperatures between +10°C to +35°C (see also chapter 2. “TRANSPORT AND STORAGE”.**

If the fume cupboard has been installed prior to building work has dried out, heating and dehumidification should be applied slowly over a number of days. Dehumidifiers should be used in preference of heaters.

If heating only will be applied, good ventilation is essential. If heating will be not suitable (summer time) provide copious natural aeration.

Where fume cupboard is installing, all surfaces, including those to be concealed by the furniture, should be decorated before setting up. Only finishing and retouches will be carried out afterward.

It is suggested to complete floor and decorated walls coverings as far as possible before fume hood is installed. Where other furniture has already been installed, it should be protected.

For further information, see the following:

- **CEN/TS 14175-5** Fume cupboards – Recommendations for installation and maintenance
- **EN 14175-4** Fume cupboards – On-site test methods
- **EN 14056** Laboratory furniture - Recommendations for design and installation
- **EN 13779** Ventilation for non-residential buildings - Performance requirements for ventilation and room-conditioning systems

## 5.4 HYDRAULIC FITTING / GASES CONNECTION

Before starting the fume hood, check the followings::


**Feeding lines** of the single fluid should be partitioned and their identification corresponding to the relevant knob color-code, according to EN 13792:2002;

**Fluid inlet pressure and delivery** are to be complying to fluid type (specially for pure gases, combustible and explosive gases);

**Drain piping** should be resistant material among poured liquids and fitted with trap;

**Special drain piping** (toxic and/or biohazard liquids) should be separated from sewer collector.

All water, gas and fluid fittings and accessories assembled on LITE fume hood manufactured according to European standards.

The remote control is provided by 2 female couplings having 3/8" gas thread; remote control knob for combustible gas is equipped with safety mechanism to avoid unintentional opening, having safety release control (open/close) and double ball valve.

The output tap is delivered with a 3/8" gas thread tang (male). The tapered hose nozzle supply according to standard.

The connection between the external control knob and the internal output tap is fulfilled by using copper tubing and mechanical seal.

The access to feeding lines is possible by removing the side cover panels, which must be carried out by skilled personnel ONLY or by technicians. Any unauthorized pilferage entails the immediate and irrevocable decadence of warranty terms and manufacturer's responsibility.

### PRECAUTIONS WITH GASES AND FLUIDS


**Avoid the use of free-flame near fuel fittings.**

Care must be taken when using sockets or flammable / explosive / toxic / corrosive gases, especially relating to the periodical seals tightness check.

The installation of alarm devices to detect gas leakages within the laboratory room is suggested.

For further informations see the chapter 10 – "SUPPLYING FIXTURE".

## 5.5 ELECTRIC CONNECTIONS

Before switching on and starting the fume cupboard, check:

**Fixing connection** to net wiring section must be complying with max input power of the assembled sockets;


**Patent voltage** on apparatus label must be the same of network supply;

**Voltage and frequency** of equipment and exhaust fan must be equivalent to network;

**Fan breaker calibration** must be complying with voltage and motor input power.


**Earthing clamp**, positioned under the working top on the back frame (see illustration) must be secured to equipotential net;

The wiring section will be conformed to hood total power electric capacity, as per label data. In case of original setting modification or alteration, the wiring feature must adapted.

Any conversion unauthorized involve the end of responsibility and guarantee by manufacturer. The modification (even partial) of the original factory preset wiring must be carried out ONLY after manufacturer's authorization and by skilled personnel or authorized technicians.


Electric supply cable must be in accordance with **IEC 60227** or **IEC 60245** standard. Use screened cable, fire retardant and halogen-free. Rated Voltage: 600/1000 V.

For further information see the chapter 10 – “SUPPLYING FIXTURE”.

### 5.5.1 ELECTRIC PANELS ON FUME HOOD

The electric components assembled on fume cupboard can be different as the control system changes (DIRECT FAN or INVERTER control or DRIVE-DAMPER control). The equipment are:

ELECTRIC SOCKETS			EXHAUST PLANT		POWER SUPPLY	
Number of sockets	I MAX (single)	I MAX (total)	FAN control	FAN power	General breaker	Electric wiring
1	8 A	8 A	INVERTER	< 1,5 kW	1P+N 16 A	3 x 2,5 mm <sup>2</sup>
2	4 A	8 A	INVERTER	< 1,5 kW	1P+N 16 A	3 x 2,5 mm <sup>2</sup>
4	2 A	8 A	INVERTER	< 1,5 kW	1P+N 16 A	3 x 2,5 mm <sup>2</sup>
1	8 A	8 A	INVERTER	> 1,5 kW	3P+N 16A	5 x 2,5 mm <sup>2</sup>
2	4 A	8 A	INVERTER	> 1,5 kW	3P+N 16A	5 x 2,5 mm <sup>2</sup>
4	2 A	8 A	INVERTER	> 1,5 kW	3P+N 16A	5 x 2,5 mm <sup>2</sup>

The electric current “I max (total)” means contemporary multiple use of sockets

Electric supply type will be conformed to fan power (if it is lower or higher than 1.5 kW).


In presence on the hood electric equipment of any **interlocked sockets**, these are **NOT PROTECTED** by the main general switch but separately powered with their own cable and supply.

## 5.5.2 MOTOR FAN

The exhaust motor fan could be connected directly to electric fume hood control panel or by ancillary circuit.


The electric wiring section must be conformed to the fan electric power, in accordance with IEC 60227 or IEC 60245 standard, using screened cable, fire retardant and halogen-free.

When distance between fume hood fitted with INVERTER DRIVE and motor fan is longer than 30 meters, an inductance must connect in-line, in order to preserve electronic control performance.

When motor fan is remote placed, e.g. on building roof, it is mandatory the installation of an **ELECTRIC SAFETY BREAKER** near the motor (in sight), in order to prevent that accidental paddle wheel will start during maintenance operation.

Anyway, see carefully the wiring diagram attached to fume hood electric power box.

**Should any doubts or misunderstandings regarding the electrical components and features, please contact the manufacturer <http://www.bicasa.it/>.**



The access to powered parts (even low voltage) is prevented by suitable safety protections (closing devices which require specific mechanical tools).

**For operator safety all the metal part of hood must be connect to earthing net. The hood earthing clamp is positioned under the working top on the back frame, near junction box.**

The V.A.V. INVERTER system **up to 1.5 Kw** motor powered will be **DELTA** connected and single-phase 230V supplied (third phase is generated by INVERTER). With power **higher than 1.5 Kw** motor will be **WYE** connected and three-phase 400V supplied. Any other different clamps connection **would be dangerous** for electric motor and V.A.V. regulation system twice.

**For further details and informations see the chapter 9 – “ELECTRIC DIAGRAMS”**

## 5.6 EXHAUST FUMES DUCTING

---

The air exhaust channel should have a section or a diameter **equivalent or bigger** than fume hood duct. The piping material must be resistant to the matters produced during fume hood processes.

Besides it is suggested:

- The use of anti acid, anti-static or non-conductive materials, fitted with gore-partitions;
- Place the fan at the end of the ducting, in order to give total depression to exhausted tubing;
- Seal junctions between the single elements forming the ducting (mandatory when exhaust fan is set at the beginning of the ducting);
- Planning and designing the plant, try to reduce number of bends (specially elbow-type) to avoid narrowing;
- Try to avoid a winding route (e.g. siphon narrowing) which may cause condensation. Just in case, an inspection drain must be provided in the lowest part of the plant;
- Calculate the pressure lost when filtering canister is in-line;
- Getting room hazard condition, suitable fire-barrier air-lock should be installed.


If few fume hood will be collected to a general duct piping (centralized aspiration) the section following the first should be at least the sum of the single plugged sections.

## 5.7 IN-LINE FILTERING

---

LITE fume hood can be equipped with filtering canister placed on aspiration ducting, preferably nearest the fan. The canister should always be installed in a position where inspection and filter replacement are convenient.

While calculating the pressure loss in the plant project, consider the filter features.

**An effective and constant maintenance and cleaning of filtering gives an additional protection to operator and to environmental and a material longer life.**

## 5.8 AUTOMATIC COMPENSATING DEVICES / MULTIPLE FUME HOODS

---

It is possible to control and to manage also articulated extraction plants, for example with a single exhaust fan servicing multiple fume hoods.

Upon request, mixed regulating systems are available, with logic control of the single fume hood according to its position or distance from exhaust fan and to the contemporary operatively factor. The automatic electronic devices can be interfaced with the building HVAC to optimize either the equipment operation, than the relevant energy saving.

## 5.9 FALSE CEILING AND DIFFUSERS


When LITE fume hood is installed in room having false ceiling or closet compartment at the top, it is necessary that the **measuring probe is connected to the ambient** in which fume hood operates.

The correct installation which allows the barometric pressure counter balance and the correct operation of the measuring system is shown on the side sketch.

Good aspiration plant and ducting are the main condition for efficient fume hood operation. When ducting crosses department or rooms partitioned in compliance with fire-prevention rule, a suitable fire-barrier air-lock should be installed. Please, consult the nearest Fire Brigade Headquarter.

The probe should be in contact with working environment in order to ensure the correct operability of measuring / check / alarm system (thus the same barometric pressure of the room from which air comes into the fume hood).

When the fume hood is fitted with filler panel up to the ceiling or installed in room having false-ceiling partitions which divide air-inlet and probe, it will be necessary the installation of a compensation grid which leads both ambient (for example a tight-tubing between the room and the probe head according to the side sketch).

## 5.10 WARNING AND PERIODIC CONTROL

Before electric box and panel opening, be sure that electric supply is **TURNED OFF**. Lack of compliance in this matter will generate very dangerous events to operators and to environmental (e.g. electrocute or fire).


After shutting off, wait few minutes to access to electronic board and circuits, looking at LED lights is power off.

**Turn ON main switch ONLY AFTER the cover and the access to electric are locked.**

When the electric supply was turned off, wait always at least 1 minute before restart.

Building air conditioned and thermal flow diffusion spots, positioned around the fume hood, near air velocity probe (it means nearer than 2 meters) may generate an abnormal system working. It get worse when room HVAC plant process defaults (unbalanced overpressure or depression condition).

The air velocity probe must be cleared and dust-free in any setting. Take away the dirty **ONLY** when the fume hood exhaust system is OFF, using compressed air (max 200Pa) or by a soft brush.

Look for electronic boxes and components abnormal color changing or smell given forth.


**DO NOT MODIFY ANY PART OR COMPONENT OF THE HOOD.**


## 6 OPERATING

### 6.1 START-UP

Turn **ON** the **GENERAL CIRCUIT BREAKER** (see 5.5.1 - 1)

Control that **THERMAL SWITCH** get **ON** (see 5.5.1 - 1)

Now LITE fume cupboard is ready to run.

### 6.2 VARIABLE AIR VOLUME SUCTION CONTROL


The fume cupboard V.A.V. control device includes master board, interface frontal board, sash opening control (optional) and air velocity probe.

The instrument is specifically designed for **AUTOMATIC EXHAUST AIR FLOW CONTROL AND SETTING**.

The setting air velocity alteration because of sash opening is automatically corrected by the electronic apparatus in order to ensure the requested performance.

The air flow delivery rate has been carried out by actuator system equipped with electric **INVERTER AC DRIVE** (for motor fan basic speed control) or with electrical **SHUT-OFF DAMPER** (on exhaust piping in-line installed).

An hot-wire sensor (temperature compensate) drives the actuator system placed in-line with frontal hood air flow and it means as anemometric device for measuring air speed.

The air speed value will be transmitted to electronic instrument that processes data input and modulates an electric signal to actuator system for keeping any set correction.

Fume hood can operate in three different ways:

- **AUTO** - with variable speed suction fan (air volume rate controlled by INVERTER DRIVE);
- **AUTO** - with fixed speed suction fan (air volume controlled rate by DAMPER);
- **MAN** - with fixed speed suction fan (air volume rate not controlled).

The instrument pilots the inverter and / or the damper by a 0/10V signal output. Damper can be controlled by 2 limit switch (open and closed damper).

**Instrument performances:**

- Face velocity setting
- Alarm parameters setup
- Buzzer setting
- Lower voltage control
- Internal light switch
- Sash auto / manual opening drive (optional)
- Sash photocell control (optional)
- Operator volumetric sensor control (optional)
- Combustible gas shut-off control (optional)
- Electric sockets shut-off control (optional)
- Inside chamber temperature detection (optional)


**6.2.1 V.A.V. SUPPLY AND INSTRUMENT LINKS****Electric supply :**

- 1x 230VAC  $\pm 10\%$ , 50 Hz – 5x20 fuse 1 Amp - INVERTER up to 1,5 kW or DAMPER;
- 3x400VAC  $\pm 10\%$ , 50 Hz – 5x20 fuse 1 Amp - INVERTER from 1,5 kW to 4,0 kW.

**Instrument links :**

- Master board – frontal panel RJ11 connect
- Master board – velocity probe RJ11 connect
- INVERTER / DAMPER pilot output
- STARTER INVERTER output
- PLC output
- Motor fan ON/OFF relay
- Lighting ON/OFF
- Alarm signal relay output
- Electric sockets relay output (optional)
- Combustible gas valve relay output (optional)
- RS232 output (remote control and board programming)
- RS485 output (sash motor drive – optional)
- Volumetric sensor input (optional)
- Sash photocell input (optional)
- Temperature probe input

## 6.2.2 V.A.V. APPARATUS

## 6.2.3 OPERATING

The control system of the LITE fume cupboard includes a master board circuit for operating performance of hood and a frontal panel with control keys and alarm device as display interface. The connection between master board and frontal panel gets on by RJ11 link.


The frontal panel is 5V supplied.

**MOTOR FAN STARTING:** to get ON or OFF the fume hood air suction system, pushing the light key **FAN** on the frontal panel.

There are two different working ways:

- **NORMAL MODE :** standard operating use of fume cupboard, with the whole available resources access.
- **SETUP MODE:** calibration of fume cupboard control, applied to set the air suction system (V.A.V.) as per operator request.

When starting operating, it checks for a VALID SETUP setting. If none would found, stop the process, restart and do INSTRUMENT CALIBRATION setup (see 6.3).

## 6.2.4 NORMAL MODE

This is the standard operating use of fume cupboard, with the whole available resources access by the FRONTAL PANEL, including:

- **Control and adjustment for base parameter**
- **Operating parameter display**
- **Operating parameter setting**
- **Alarm management**
- **Optional performance and control**

## 6.2.5 BASE PARAMETER

### LIGHTING

The LIGHT key switches ON or OFF the ceiling fluorescent lamp.

- YELLOW LED indicates the instrument is supplied and light is OFF.
- GREEN LED indicates the light is TURNED ON.

### AIR SUCTION

The FAN key switches ON or OFF the motor fan connected to the hood.

- YELLOW LED indicates the instrument is supplied and fan is OFF.
- GREEN LED indicates the fan is TURNED ON.

### AIR VELOCITY DISPLAY

The display normally shows the frontal opening air velocity with 3 figures (0,00 meter per second).

There are three different air suction capacity operating processes, selectable by specific duties:

- **Automatic air suction mode**
- **Manual air suction mode**
- **Emergency air suction mode**

### AUTOMATIC AIR SUCTION MODE

The frontal opening air velocity has been AUTOMATICALLY steady maintained by the instrument control, by the anemometric probe and by the power control system (INVERTER or DAMPER), even if sash opening has been changed. The automatic process is typically in ordinary working conditions.

- GREEN LED signal means the AUTOMATIC mode is TURNED ON.

The frontal air rate must be set by the “FACE VELOCITY” program (see 6.2.6).

### MANUAL AIR SUCTION MODE

The frontal opening air velocity has been MANUALLY set by operator. The manual mode replaces the automatic function: the air flow rate is steady maintained but the frontal air velocity changes if sash opening will be modified. The manual process is normally used in specific working conditions, when it is mandatory to operate with a particular or different frontal air velocities.

● GREEN LED signal means the MANUAL mode is TURNED ON.

The flow rate shows on the figures display; it must be set by the keys  and .

### EMERGENCY AIR SUCTION MODE

Pushing the RED key “EMERGENCY” during process, the frontal opening air velocity will turned up to MAXIMUM rate, regardless the just setting condition. Pushing a second time the key, the operating mode will go back to previous setting. The emergency process is used in dangerous or safety risk conditions, when it is advisable to run the maximum fan exhaust capacity.

● RED LED signal and acoustic BEEPER indicate EMERGENCY mode is TURNED ON.

To stop the emergency state and return to previous condition, push again the key till the RED LED is OFF.

## 6.2.6 WORKING PARAMETERS

The hood electronic device allows to set the single control parameters and linked alarms. The available programming data are the following:

- FACE VELOCITY
- SASH TIMER
- ALARM VELOCITY
- ALARM TEMP

**PUSH the “MENU” key** to see on the DISPLAY the related data, switching on the own LED signal. The lighting side ● GREEN LED will indicate the programming function. Push time after time the key to access to any parameter.

The available SETTING data are:

- **FACE VELOCITY** Air SUCTION SPEED value on the frontal sash opening, benchmark of the hood automatic set-point in use. The rate may be fitted in the range from 0.10 to 0.90 m/s. The default set point is 0.50 m/s.

- **SASH TIMER** (OPTIONAL) Delayed TIME of the automatic motor driven sash closing. This function is available only when the fume cupboard is fitted with motor driven sash and with operator presence sensor. The rate may be fitted on 0 / 2 / 4 / 6 minutes (0 = disable).
  
- **ALARM VELOCITY** SUCTION SPEED MINIMUM ALARM value on the frontal sash opening. The rate may be fitted in the range from 0.20 to 0.80 m/s. The default set point is 0.30 m/s.
  
- **ALARM TEMP** (OPTIONAL) MAX TEMPERATURE value of the internal chamber air flow. This function is available only when the fume cupboard is fitted with optional temperature probe. The rate may be fitted in the range from 10°C to 60°C. The default set point is 30°C.

## 6.2.7 WORKING PARAMETERS SETTING

Operator may set the available working data following standard procedure.

	“SEL” KEY to start the programming mode
	“MENU” KEY to choose the changeable data. Push time after time the key to scroll any parameter. The selected data will be pointed by the <span style="color: green;">●</span> GREEN LED lighting on.
 	“PLUS” and “MINUS” KEY to adjust the displayed rate.
	“MENU” KEY to validate the new data entering.
	“SEL” KEY to get out from the programming mode .

## 6.2.8 ALARM MANAGEMENT

The safety operating of the fume cupboard has been available by electronic control facilities apparatus. The alarm signal have been ACOUSTIC and LIGHTING device.

	<p>The buzzer ACOUSTIC signal is installed on master board circuit and it could be silenced (for 2 minutes timing) pushing at the same time “PLUS” and “MINUS” keys. Upon expired the period, if the failure condition was still present, the buzzer will ring again.</p> <p>The LIGHTING signal will stay even ON up to alarm failure will get down.</p>
---	---

The available alarm devices are the following:

 <b>ALARM VELOCITY</b>	<p>● RED LED signal and acoustic BEEPER indicate the malfunction.</p> <p>It had been when the air flow velocity failed compared to operating set value, for a period of 30 seconds or more.</p>
 <b>ALARM TEMP</b>	<p>● RED LED signal and acoustic BEEPER indicate the malfunction.</p> <p>It had been when the internal chamber temperature was higher than the pre setting point. The utility is available when optional kit is fitted.</p>
 <b>ALARM SASH UP</b>	<p>● RED LED signal and acoustic BEEPER indicate the malfunction.</p> <p>It had been when the sash opening goes override, beyond the 50 cm of height. The utility is available when optional kit is fitted.</p>

**Alarm condition event is ALWAYS a potential RISK CONDITION for operator. Carefully detect as soon as possible the generating wreck and do as anything as possible to remove.**

If necessary, shut off hood and call TECHNICAL SUPPORT (see the chapter 1.1)

### LOW-VOLTAGE TROUBLE

In case of electric black out, the running function will be protected and saved until supply will come back. Including the EMERGENCY condition.

## 6.2.9 OPTIONAL PERFORMANCES

The hood electronic device allows to fit several OPTIONAL devices.

- **SASH MOTOR DRIVE**
- **COMBUSTIBLE GAS VALVE CONTROL**
- **ELECTRIC SOCKETS SUPPLY CONTROL**
- **“SASH UP” ALARM DETECTOR**

### **SASH MOTOR DRIVE:**

by the optional supplementary circuit is possible to control the electrical opening and closing of the sash, pushing the “SASH”  UP and  DOWN keys. The device fitted with motor drive, transmission clutch, photocell safety barrier and operator presence sensor (delayed closing).

- YELLOW LED color will indicate the utility carry on.
- GREEN LED color means the device is working.

### **COMBUSTIBLE GAS VALVE CONTROL:**

by the optional supplementary circuit is possible to control the combustible gas inlet by electric valve control, just when suction fan runs. The “GAS” key will manually open (only with air suction) or close the valve.

- YELLOW LED color will indicate the utility carry on.
- GREEN LED color means the device is working.

### **ELECTRIC SOCKETS SUPPLY CONTROL:**

by the optional supplementary circuit is possible to control the electric sockets supply control, just when suction fan runs. The “FLASH” key will manually switch ON (only with air suction) or switch OFF the supply.

- YELLOW LED color will indicate the utility carry on.
- GREEN LED color means the device is working.

### **“SASH STOP DEVICE” ALARM DETECTOR:**

According to EN 14175 the fume cupboard is fitted with a sash stop device, limiting height in operational sash opening.

If needed, it is available an optional ACOUSTIC and LIGHTING alarm device to indicate the opening is overridden.


## 6.3 INSTRUMENT CALIBRATION


This section of the manual is reserved to instrumentation **SKILLED** specialist and technician.

The wrong use of calibration set could do abnormal functioning of the fume cupboard and done **HAZARDOUS** condition.

On this facility mode it is possible to calibrate the control system of the hood according to environmental and

The calibration **MUST** be done at the end of installation process and any time the anemometric probe has been replaced, to guarantee the right performance of the hood.

Setting of MENU parameters designed for highly qualified personnel. An accidental or not controlled change of values may be dangerous for fume cupboard performances and operator safety.

### 6.3.1 CALIBRATION MENU

Before starting the calibration menu control that the “MENU”  key led signal be OFF.

Push at the same time “PLUS”  and “SASH UP”  keys.

On DISPLAY will appear the word “Cod”.

Use the “PLUS”  and “MINUS”  keys to input the password (default value is **010**).

Push “LIGHT”  key to confirm data.

At this step on DISPLAY will appear the following function entries:

“OFF”	To adjust the frontal air velocity offset
“UEL”	To set the inverter minimum operating frequency
“CAL”	To calibrate the frontal air velocity measurement

Use the “UP”  and “DOWN”  keys to scroll the list.

## 6.3.2 “CAL” - FRONTAL AIR VELOCITY MEASUREMENT CALIBRATION

The frontal air velocity measurement calibration must be done with 2 benchmarks (LOW and HIGH).

Raise the sash **up to about 10 cm of height** from air-foil or from countertop level. Hold the anemometric or similar tool measuring probe on steady support, approximately in the middle of sash opening suction front, approximately at 30 cm outside.


Avoid turbulence or interfering flow (close doors and windows) and operator obstruction in front of sash opening.

Upon accessing to the CALIBRATION MENU, pushing UP  and DOWN  keys select the “CAL” function; then push LIGHT  key to run.

On the display it will appear:

- **001** - LOW CALIBRATION
- or **002** HIGH CALIBRATION

Use UP  and DOWN  keys to select the mode and push LIGHT  key to run.


### LOW CALIBRATION (001)

On the display it will appear a LOW air velocity value (e.g. 0.30). With UP  and DOWN  keys do change of fan electric motor speed controlled by the INVERTER (the variation could request few seconds to go) until the benchmark corresponding to the outside anemometer rate is fitted, or rather a frontal velocity equal to **0,30 m/s**.

Wait for about 2 minutes. Do a fine calibration set using the “PLUS”  and “MINUS”  keys, compared to the external anemometer stabilized rate.

Upon the system settled, push LIGHT  key to confirm data. Wait for few moments for the data processing up to the previous MENU reset.

Go to HIGH CALIBRATION 002 using UP  and DOWN  keys, then push LIGHT  to run.

### HIGH CALIBRATION (002)

On the display it will appear a HIGH air velocity value (e.g. 0.70). With UP  and DOWN  keys do change of fan electric motor speed controlled by the INVERTER (the variation could request few seconds to go) until the benchmark corresponding to the outside anemometer rate is fitted, or rather a frontal velocity equal to **0,70 m/s**.

Wait for about 2 minutes. Do a fine calibration set using the “PLUS”  and “MINUS”  keys, compared to the external anemometer stabilized rate.

Upon the system settled, push LIGHT  key to confirm data. Wait for few moments for the data processing up to the previous MENU reset.

At the end of process, push at the same time “PLUS”  and “SASH UP”  keys to return to calibration main MENU.

### 6.3.1 “OFF” – AIR VELOCITY MEASUREMENT COMPENSATION

After choosing the “OFF” mode on the main CALIBRATION MENU, push at the same time the UP  and DOWN  keys and LIGHT  to run. On the display it will appear the CORRECTION READING VALUE (e.g. .00).

Using the “PLUS”  and “MINUS”  keys to select the OFFSET value and then push LIGHT  key to confirm and to return to main MENU.

After data adjustment, the whole air velocity values processed by the instrument will be changed at the same number (more or less) to the actual one. It means that, if the instrument display shown previously 0.47 and the set “OFF” value has been set at .02, the new figure on display will be 0.49; if the “OFF” value has been set at - .02 the new figure on display will be 0.45. The rectification is standing up to new configuration.

**ADJUST / LINEUP THE DISPLAY VALUE ONLY IF IT IS REALLY NECESSARY.**

At the end of process, push at the same time “PLUS”  and “SASH UP”  keys to return to calibration main MENU.


### 6.3.2 “UEL” – INVERTER MINIMUM FREQUENCY RATE

This function set the inverter minimum frequency value share, that is the rpm fan motor minimum speed, fixed by manufacturer. The rate could be solely INCREASED up to 100% (it could not be less then manufacturer base).

Choose the “UEL” mode on the main CALIBRATION MENU, push at the same time the UP  and DOWN  keys and LIGHT  to run. On the display it will appear the default MINIMUM FREQUENCY share (e.g. 20).

To increase and refine the value, use the “PLUS”  and “MINUS”  keys,

Push LIGHT  key to confirm and to return to main MENU.

**TO BACK OUT FROM CALIBRATION MENU AND RESTORE NORMAL OPERATING OF THE HOOD, SHUT OFF THE SYSTEM BY MAIN GENERAL SWITCH AND RESTART.**

## 6.4 AIR VELOCITY MEASURING PROBE

The air velocity gauge is done by hot-wire anemometric probe. Hot wire anemometers use a very fine wire (on the order of several microns) electrically heated up to some temperature above the ambient. Air flowing past the wire has a cooling effect on the wire. As the electrical resistance of most metals is dependent upon the temperature of the metal, a relationship can be obtained between the resistance of the wire and the flow speed.

Probe is linked to master board circuit by RJ11 cable; power supply. 5V.

Presence of anemostat or air conditioned diffuser close to the hood ceiling and the probe (less than 2 meters) can do abnormal functioning of control system, as well as the barometric room pressure balancing failure.

## 6.5 INVERTER SSD DRIVES - 650 SERIES

The 650 series is a family of AC drives that provide a no-fuss, open-loop AC motor control applications. The 650 has single key selectable pre-programmed applications so set-up is quick and easy without superfluous complications. Powers are available from 0,25KW up to 7.5KW on 230V single-phase or 400V three-phase.

The 650 series inverter starting and operating procedure had been designed for the easily parameter configuration and calibration.

Configuration and parameter setting carried out in warehouse.

INVERTER SET CONFIGURATION		
Parameter	Value	SSD 650 Settle
P1	1	application
P2	50Hz	frequency max
P3	10%	minimum speed (related to cm 180 hood)
P4	5	acceleration rate (seconds)
P5	5	deceleration rate (seconds)
P6	2 ÷ 5	see motor to drive (4÷)
P7	50Hz	normal frequency
P8	10,0	speed (when "joy" is high)
P9	0	stopping mode
P11	1	identifying the charge type
P12	0	current limit (150% for 30 seconds)
P13	5%	low frequency power flow

INVERTER PARAMETER CHANGING		
1.	Push for any seconds the <b>M key</b>	It will appear <b>Djoy</b>
2.	Push <b>▲</b>	It will appear <b>PAR</b>
3.	Push the <b>M key</b>	It will appear <b>P1 ... P13</b> depending on last access
4.	Select the parameter <b>PX</b> with <b>▲ ▼</b>	
5.	Push again the <b>M key</b>	It will appear the related value
6.	Set the parameter value with <b>▲ ▼</b>	Confirm the value with <b>E key</b>
7.	Set the following with <b>▲ ▼</b>	Repeat the sequence from <b>point 4.</b>
At the end of setting, push any time the <b>E key</b> up to <b>READY</b>		

## 6.6 CIRCULAR SHUT-OFF DAMPERS – ELECTRONIC ACTUATOR

The shut-off dampers are employed to modulate volume flows in ventilating and air conditioning systems by V.A.V. instrument control. The control then let the system know to reduce or increase the volume of air that needs to be exhausted.

Damper is PVC tubing made and adjustment can be carried out by a 24VAC electric gear motor.

## 6.7 FUMES SUCTION AND EXHAUST

A good aspiration plant and ducting, are the main conditions for the chemical fume hood operation. Specifically, the fume hood connection to a system having tailored processing parameters, is a milestone, since the air-volume changes and air-face velocity are strictly connected to exhaust plant setting.

LITE fume hoods – BENCH TOP model (with countertop at cm 90 of height)							
FACE VELOCITY	0,25 m/s		0,50 m/s		0,65 m/s		ΔP
	Output air	Volume per min	Output air	Volume per min	Output air	Volume per min	
LITE 1.2	470 m <sup>3</sup> /h	8	1030 m <sup>3</sup> /h	16	1210 m <sup>3</sup> /h	21	100 Pa
LITE 1.5	590 m <sup>3</sup> /h	8	1300 m <sup>3</sup> /h	16	1530 m <sup>3</sup> /h	21	112 Pa
LITE 1.8	710 m <sup>3</sup> /h	8	1580 m <sup>3</sup> /h	16	1840 m <sup>3</sup> /h	21	125 Pa
LITE 2.4	940 m <sup>3</sup> /h	8	2100 m <sup>3</sup> /h	16	2450 m <sup>3</sup> /h	21	125 Pa

Reported values are referred to max sash opening of cm 50 with tolerance set forth by EN 14175 standard.

With V.A.V. system control, when sash is closed, values are drastically reduced.

If the fume hood is equipped with acoustic and visual alarm device, eventual failures will be detected and a buzzer will be activated. To stop the alarm, push relevant key on instrument panel.

**It is strictly forbidden using the fume hood without ventilation thus without mechanical convection.**

**Air shortage within the room means a decrease of fume hood performances and less safety for the operator.**

## 6.8 SASH


The frontal sash is employed to easy access to working hood chamber. It consists of a vertical moving frame including horizontal sliding glasses.

A safety mechanical lock according to EN 14175 standard stops sash opening at 50 cm of height. This limit should be override only during maintenance, when introducing large instrument within working chamber or for specific reasons.

When the sash come back lower, the safety lock will be automatically restored without the need of mechanical tools.

Sliding glasses opening do not have safety lock.

During chemical process **DO NOT OPEN** the vertical sash and the sliding glasses **AT THE SAME TIME**.

Safety lock must be removed releasing the levers positioned on the right side under the lifting handle **1**.

**Keep sash closed whenever possible.**

**Avoid to keep twice glasses on the same sash frame side.**

## 7 CARE AND MAINTENANCE

### 7.1 GENERIC ADVICES

---

A correct pollution guarantees the operator safety and the apparatus shelf-life. The fume cupboard outside parts may be treated with usual cleaning agent, as well as the internal chamber when in the present of inorganic waste, weak acid spot, chalky residual.

With organic reagents use specific disinfectant cleaning matter, carried out works suitable.

In case of very hazardous pollution risk, chemical or biological, utilize sterilize or decontaminator product **as long as chlorine-free and/or chlorine-derivate-free.**

### 7.2 WORKING AREA ACCESS

---

**According to EN 14175 - Part 2: Safety and performance requirements - 7.3.1 and 7.3.2 paragraph,** the operational sash incorporates a stop device to prevent the sash being opened more than the maximum operational sash opening (500 mm).

It shall not be possible to override the stop without a deliberate act on the part of the operator. The design of the sash stop shall ensure that the stop is automatically reset when the sash is returned to a position less than the maximum operational opening.

- It is better to override the sash stop limit only for maintenance or to place large volume apparatus inside the working chamber: anyway **do not override sash stop during analysis or chemical test performances nor with exhaust fan shut off.**
- It is appropriate **to keep sash closed whenever possible:** the frontal air velocity changes proportionally by the sash opening. When the sash is higher opening, the lower is the air frontal velocity.

### 7.3 SASH AND COUNTERWEIGHT MAINTENANCE

---

Test periodically the vertical sash device scrolling up and down, checking the stainless steel ropes status and the counterweight balancing. In case of ropes repairs as well as in case of introducing large volume apparatus and equipment, override the sash stop (totally opening) unblocking the control device and activating the 2 small lever placed under the sliding-glass handle

**Do not maintenance operation and/or do not override sash stop during analysis or chemical test performances nor with exhaust fan shut off.**

Keep the sash as high as possible as the counterweight get to end point, on resting position.

The ropes are simply accessible:

- In the front part (behind the sash). On the internal face of lateral sash frame catching ropes, their anchor bolts and their sliding regulating skid. The plastic skates couple with metal plates is setting by adjustable buttonholes with lateral fixing. Making the modification, be sure that both vertical sash free-running.
- In the upper part (ceiling). Check the pulleys status and the ropes driving into. The pulleys must be free revolving, the cord must be precise housed in the gorge of the pulley, inside the bracket cover.
- In the rear part (internal baffle). The terminal parts of counterweight are accessible removing the 2 bottom inside panels (screwed). It is possible to check the correct hook up, to couple or uncouple the rope, adjust the rope stretching by the fixing screw.  
Le funi sono raggiungibili con facilità:

**The maintenance operations must be performed by skilled personnel.**

## 7.4 TO AVOID

---

- The continuous powder and/or abrasive sponge use may degrade the polished and painted surfaces.
- The hydrofluoric acid and derived from, used to clean, corrupt glass and polished surfaces.
- The hydrochloric acid and derived from, used to clean, corrupt stainless steel surfaces.
- The free flames may damage to paint and plastic laminated colours.
- The total opening sash working is dangerous: keep it closed whenever possible, to minimize operator risk.

**It is very important that these processes should be carry-out by skilled personnel.**


## 8 SPARE PARTS

PRODUCT	CODE	HOOD TYPE
<b>GLASS</b>		
Side shoulder glass panel	FAL-220330	LITE 1.2 / 1.5 / 1.8
Sliding glass door	FAL-220332	LITE 1.2
Sliding glass door	FAL-220333	LITE 1.5
Sliding glass door	FAL-220334	LITE 1.8
Illuminating lamp glass shield	FAL-190603	LITE 1.2 / 1.5
Illuminating lamp glass shield	FAL-190605	LITE 1.8
<b>LIGHTING</b>		
Fluorescent tube - D26 white – 18W	FAL-120143	LITE 1.2 / 1.5
Fluorescent tube - D26 white – 36W	FAL-120144	LITE 1.8
<b>SASH</b>		
Counterweight skid	FAL-180272	
Sash skid	FAL-160845	
Rope stretching	FAL-200357	
Revolving pulley	FAL-030064	
Sash stop device spring	FAL-130194	
<b>V. A. V. " CONTROLLER .</b>		
Frontal panel	FAL-191511	
Master board	FAL-191505	
Hot wire sensor (temperature compensate)	FAL-191510	
RJ11 junction board	FAL-191512	
RS232 Extension kit	FAL-110173	
Frontal panel to master board connection kit	FAL-110174	
Single phase inverter (220-240V) 0,75 kW	FAL-090021	
EMC outside filter - Class B 0,75 kW	FAL-060296	
Single phase inverter (220-240V) 1,5 kW	FAL-090012	
EMC outside filter - Class B 1,5 kW	FAL-060297	


## 9 ELECTRIC DIAGRAMS


### 9.1 AIR SUCTION CONTROL DIAGRAM – POWER UP TO 0,75 kW


## 9.2 AIR SUCTION CONTROL DIAGRAM – POWER UP TO 1,5 kW


## 9.3 GENERAL BLOCK DIAGRAM

**10 SUPPLYING FIXTURE**

1800	250
1500	250
1200	200
<b>L</b>	<b><math>\phi</math></b>

**ON WALL**

## FITTINGS

- S** = WASTE –  $\phi$  50 / 60 mm acid-proof pipe
- A** = WATER –  $\phi$  3/4" or 1/2" zinc-plated steel pipe (according to number of taps) provided with main valve
- DEIONIZED WATER –  $\phi$  1/2" pipe made of stainless-steel or PVC
- G** = GASES –  $\phi$  1/2" zinc-plated steel or  $\phi$  16 mm copper pipe provided with main valve
- E** = ELECTRICAL – water-tight terminal box for single and / or three-phase + N power line with earth lead (according to services)

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## 11.4 ANSI/ASHRAE 110-1995 CERTIFICATIONS

ANSI/ASHRAE 110-1995 Flow Containment & Instrumentation

**CERTIFICATE OF TYPE TESTING  
IN ACCORDANCE WITH ASHRAE 110-1995**

CERTIFICATE & REPORT NO: INV/ASHRAE110/411  
DATE: 19th November 2007

**Fume Cupboard Manufacturer:**  
Bicasa spa  
I-20044 Bernareggio (Milano)  
Viale delle Industrie 33  
ITALY

**Fume Cupboard Type:**  
LITE, 1.2m wide bench-type combination sash

Vertical test opening width:	1110mm
Vertical test opening height:	500mm
Horizontal test opening width:	480mm
Horizontal test opening height:	840mm
Average face velocity at vertical sash opening:	0.50m/s

This is to certify that the fume hood described above has been type-tested in accordance with ASHRAE 110-1995 and resulted in velocity and containment performance given in the corresponding test report.

Certified by: Dr A F Bicen  
Invent UK Ltd, 85 Southdown Road, Harpenden, Herts AL5 1PR  
Tel: 01582 461666 : Fax: 01582 460077 : E-mail: info@invent-uk.com

ANSI/ASHRAE 110-1995 Flow Containment & Instrumentation

**CERTIFICATE OF TYPE TESTING  
IN ACCORDANCE WITH ASHRAE 110-1995**

CERTIFICATE & REPORT NO: INV/ASHRAE110/412  
DATE: 19th November 2007

**Fume Cupboard Manufacturer:**  
Bicasa spa  
I-20044 Bernareggio (Milano)  
Viale delle Industrie 33  
ITALY

**Fume Cupboard Type:**  
LITE, 1.5m wide bench-type combination sash

Vertical test opening width:	1410mm
Vertical test opening height:	500mm
Horizontal test opening width:	500mm
Horizontal test opening height:	840mm
Average face velocity at vertical sash opening:	0.50m/s

This is to certify that the fume hood described above has been type-tested in accordance with ASHRAE 110-1995 and resulted in velocity and containment performance given in the corresponding test report.

Certified by: Dr A F Bicen  
Invent UK Ltd, 85 Southdown Road, Harpenden, Herts AL5 1PR  
Tel: 01582 461666 : Fax: 01582 460077 : E-mail: info@invent-uk.com

ANSI/ASHRAE 110-1995 Flow Containment & Instrumentation

**CERTIFICATE OF TYPE TESTING  
IN ACCORDANCE WITH ASHRAE 110-1995**

CERTIFICATE & REPORT NO: INV/ASHRAE110/413  
DATE: 19th November 2007

**Fume Cupboard Manufacturer:**  
Bicasa spa  
I-20044 Bernareggio (Milano)  
Viale delle Industrie 33  
ITALY

**Fume Cupboard Type:**  
LITE, 1.8m wide bench-type combination sash

Vertical test opening width:	1710mm
Vertical test opening height:	500mm
Horizontal test opening width:	500mm
Horizontal test opening height:	840mm
Average face velocity at vertical sash opening:	0.50m/s

This is to certify that the fume hood described above has been type-tested in accordance with ASHRAE 110-1995 and resulted in velocity and containment performance given in the corresponding test report.

Certified by: Dr A F Bicen  
Invent UK Ltd, 85 Southdown Road, Harpenden, Herts AL5 1PR  
Tel: 01582 461666 : Fax: 01582 460077 : E-mail: info@invent-uk.com

**COPYRIGHT**

All the technical information content into this manual, as per drawings, figures, picture and so long, are protected by copyright international law and cannot be used, copied, multiplied, send or given way to anybody without written authorization of the manufacturer.

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BIOLOGICAL SAFETY CABINET

# APPLICATIONS BIO II A

BIO II A cabinet provides personnel, environmental and product protection in applications involving biological agents. These cabinets have been designed according to EN 12469 and are meeting or exceeding requirements of NSF 49, BS 5726, NFX 44-201 and DIN 12950.

Each cabinet is factory-certified by TÜV-GS registered laboratory to meet the specified performance requirements.

## CONTROL SYSTEM

Main switch with key and microprocessor control system. Alphanumeric luminescent display providing following data:

- Exhaust air flow  $m^3/h$ .
- Laminar flow air velocity  $m/s$ .
- Cabinet temperature  $^{\circ}C$ .
- Fans elapsed hour meter.
- U.V. lamp elapsed hour meter.
- HEPA filters last change date.

Safety and regulation systems controlled by the microprocessor. Visual and auditory alarms shown on the display:

- Low exhaust flow.
- Low downflow air velocity.
- Exhaust fan malfunction.
- Front window opened.

Control system is also providing safety features as:


- U.V. interlocking when front window is open.
- Air flow and lighting interlocking when U.V. lamp is on.
- Impulsion and exhaust fans providing a double safety system: in case of eventual stop of the exhaust fan, the impulsion

one is automatically stopped. In case of eventual stop of the impulsion fan, the exhaust one increases its speed in order to maintain underpressure conditions in the working space.

-Microprocessor compensates for gradual increase in filters resistance: inflow and downflow are maintained in any condition.

# GENERAL FEATURES

BIO II A cabinet is 70% recirculating down flow enclosure Class II A2 equipped with two high efficiency HEPA filters:

- HEPA filtered laminar airflow Class 10, which is recirculated in the work zone, creates an ultra clean work environment for product protection. A room air barrier at the front opening prevents escape of pathogens into the worker's area.

- Exhaust HEPA filter provides safety to environment.

Stove-enamel coated steel cabinet. Easy access for working and cleaning. Working space in stainless steel AISI 304. Segmented working surface. Laminated safety front window with gas springs. Two HEPA/ULPA filter stages H 14 efficiency 99.999% for particles  $>0.3 \mu$  easily accessible for replacement.

Vertical down flow filter, exhaust filter and fans accessible from the front provide easy maintenance.

Lighting system outside working chamber makes service easier.

Possibility to connect an output relay for extend fan activation (ducted cabinets).

Special model with lead protection for radioisotopes preparation.

Front door in tempered glass with gas springs

Standard delivery includes U.V. lamp, electrical socket and two stopcocks


# TECHNICAL DATA

FEATURES	UNIT	BIO II A/P	BIO II A	BIO II A/M	BIO II A/G
External dimensions (crated)	DxWxH mm	850x970x1630	850x1420x1640	960x1700x1630	960x1990x1670
External dimensions (uncrated)	DxWxH mm	760x852x1410	760x1310x1410	760x1615x1410	760x1920x1410
Weight (crated/uncrated)	Kg	180/160	234/200	274/235	329/270
Power	kW/A	1.2	1.8	1.8	2.9
Electrical supply	V/Hz	230/50 single ph (other supplies available)			
Chamber dimensions (DxWxH)	mm	580x732x700	580x1190x700	580x1495x700	580x1800x700
Laminar downflow velocity	m/s	0.40	0.40	0.35	0.40
Downflow	m <sup>3</sup> /h	605	995	1100	1500
Exhaust flow	m <sup>3</sup> /h	390 (Ø 160mm)	430 (Ø 160mm)	550 (Ø 200mm)	770 (Ø 200mm)
Inflow velocity	m/s	0.72	0.54	0.51	0.6
Light intensity	lux	≥ 760	≥ 930	≥ 1100	≥ 1100
Noise level	dBA	< 58			
Filters HEPA/ULPA H14	%	Efficiency > 99.999% (class 10)			


Standard configuration: BIO II A work station including UV-light, 2 stopcocks (one with safety solenoid valve), electrical socket in working area and electrical cable with plug (specify type when ordering).

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**CASELLA ESPAÑA S.A.**